



SUCTION-/IRRIGATION SET; DISPOSABLE



SPECIAL SALES OFFER

until 31/10/2018 (while stock last)

SUCTION-/IRRIGATION SET; DISPOSABLE Ø 5.0 mm **FL** 330 mm (13")

Description	Art. No.
<ul style="list-style-type: none"> ▶ 3.0 m tubing; 2 spikes ▶ stainless steel cannula; Ø 5.0 mm; 330 mm ▶ disposable ▶ sterile ▶ 25 pcs./box 	521-090-05

This suction-/irrigation set offers an integrative system of disposable tubing and valve in combination with the reusable cannula which allows the physician to comfortably and efficiently irrigate and aspirate during and after the surgical procedure.

- ▶ Flow controlling trumpet valve
- ▶ Plunger caps are color-marked and contour-coded
- ▶ Equipped with twin irrigation spikes and blocking clamps
- ▶ Fitted with atraumatic tip with suction relief holes
- ▶ Smooth and reliable operation control for irrigation and aspiration
- ▶ Ergonomic design